

The progress of CARR's neutron imaging facilities

Hongli Wang, Songbai Han, Lijie Hao, Meimei Wu,
Dongfeng Chen, Yuntao Liu, Kai Sun

2009. 10. 28



China Advanced Research Reactor (CARR)



Key Parameters

- 60 MW
- Max undisturbed thermal neutron flux ($n \cdot \text{cm}^{-2} \cdot \text{s}^{-1}$)
 - 1×10^{15} (at reactor core)
 - 8×10^{14} (at heavy-water reflector)
- 19.75 wt% U^{235} enrichment









Some examples of neutron imaging applications


visible light sees only the








X-ray (150keV) cannot see through



Gamma-ray (1.2MeV) shows too low contrast




Neutron (25meV) reveals explosive interior details


Neutron imaging reveals the rubber o-ring and adhesive inside metallic body

Neutron Tomography



interior heat conductor

Thermal and Cold Neutron Imaging Facility Location





Thermal Neutron
Imaging Facility

Cold Neutron
Imaging Facility

Thermal Neutron Imaging Facility

Dr. Lijie HAO, 3 staffs




Design Parameters of Thermal Neutron Imaging Facility

- **CARR reactor: 60 MW**
- **Maximum of thermal neutron flux at sample : $3 \times 10^8 \text{ n cm}^{-2} \cdot \text{s}^{-1}$**
- **Source-object distance: 10 m**
- **Measurement position (L): 3m, 6m, 9m**
- **L/D: 50~9000**
- **Cd-ratio : >100**
- **n/γ : $>1.0 \times 10^6 \text{ n cm}^{-2} \cdot \text{mR}^{-1}$**
- **Beam size at sample position : 250 mm x 250 mm**
- **Desired resolution : 0.15mm,**

The parameters of Thermal Neutron Imaging Facility



D (cm)	L (cm)	L/D	Neutron flux at sample (n cm⁻²·S⁻¹)	Cd-ratio
6	900	150	2.3×10^8	>100
	600	100		
	300	50		
4	900	225	1.0×10^8	>100
	600	150		
	300	75		
2	900	450	2.5×10^7	>100
	600	300		
	300	150		
1	900	900	6.3×10^6	>100
	600	600		
	300	300		
0.1	900	9000	6.3×10^4	>100
	600	6000		
	300	3000		

Flight Tube System




Beam Limiter

The beam limiter is installed at the end of the flight tube. It consists of four boron aluminum plates that can be driven independently into the neutron beam




Sample tables

The large rotation table




The small rotation table



CCD Camera





Shielding Gate




Cold Neutron Imaging Facility


Dr. Meimei WU, 3 staffs




Diaphragm




Flight tube





flight tube

Beam limiter





Build a tight cooperation in the near future!

Thank you for your attention!

